PROTOCOL FOR ISOLATION OF PHYTOPHTHORA

(Erwin & Ribeiro, 1996)

- Wash the diseased tissue thoroughly in running tap water.
- ➤ Surface sterilize it with 0.1% mercuric chloride for 1min (time of exposure will change with the type of the tissue).
- ➤ Wash it with three changes of sterile distilled water.
- > Blot it dry on a sterile filter paper.
- ➤ Transfer it to selective medium for isolation of *Phytophthora* (PVPH medium)
- ➤ Incubate it at room temperature for 5 days.